

## PS Capture™ Exosome ELISA Kit (Anti Mouse IgG POD)

Wako

Wako

96 tests

Rating: Not Rated Yet

### Price

Sales price 6510 ?

Discount

*On Order*

★★★★★

[Ask a question about this product](#)

Manufacturer [FUJIFILM Wako](#)

### Description

- Other information
- [SDS](#)
- Product specification sheet
- [English](#)
- Package Insert
- [English](#)

This product is for research use only. Do not administer it to human.

This kit is an enzyme immunoassay reagent that can be used for the qualitative analysis of extracellular vesicles purified from cell culture supernatants and body fluid specimens and quantitative analysis of extracellular vesicles in cell culture supernatant specimens. Extracellular vesicles with optional marker protein on the surface can be detected with high sensitivity by using a solid phase immobilized protein that binds specifically to phosphatidylserine (PS) of extracellular vesicles and an arbitrary detection mouse monoclonal primary antibody.

### Features

1. High sensitivity detection of extracellular vesicles by Tim 4 affinity method for phosphatidylserine
2. Quantitative measurement of extracellular vesicles in culture supernatant (requires extracellular vesicle standard purified with MagCapture™ Exosome Isolation Kit)
3. Sensitivity of 50 to 1000 times or more of the western blot method
4. Detection control Anti-CD63 antibody attached as primary antibody

### Reviews

There are yet no reviews for this product.