

## Viral Nucleic Acid Extraction Kit II (300prep)



300T

Rating: Not Rated Yet

**Price**

Sales price 4305 ?

Discount

*On Order*



[Ask a question about this product](#)

Manufacturer [Favorgen](#)

Description

- [Specification](#)
- [Component](#)

Principle: spin column (silica membrane) Sample: 150 µl cell-free fluid such as serum, plasma, body fluid and cell cultured supernatant  
Fragment size: 100 bp ~30 kb Recovery rate: 70 ~ 90 % Binding capacity: 30 ug Elution volume: 40 ~ 50 µl Operation time: 20 min

FAVVK 001(50 preps)FAVVK 000-1(4 preps\_sample)FAVVK 001-1(100 preps)FAVVK 001-2(300 preps) VNE Buffer 35 ml1.8 ml x 270 ml200 ml Wash Buffer 1 (concentrate) 22 ml0.9 ml x 244 ml132 ml Wash Buffer 2 (concentrate) 20 ml1.5 ml20 ml x 250 ml x 2 RNase-free Water 6 ml0.5 ml12 ml20 ml Carrier RNA 0.4 mg-----0.8 mg2.2 mg VNE Column 50 pcs4 pcs100 pcs300 pcs Collection Tube 100 pcs8 pcs200 pcs600 pcs Elution Tube 50 pcs4 pcs100 pcs300 pcs User Manual 1111

**Reviews**

There are yet no reviews for this product.