

Corning DeckWorks 1 - 300uL Pipet Tips, Graduated, Reload Starter Kit, Natural, Non-Sterile, Polypropylene, 96 Tips/Deck, 10 Decks/Pack, 960 Tips/Pack, 1 Pack/Case. and 10 Empty Ra



960pcs/CS

Rating: Not Rated Yet

Price

Sales price 548 ?

Discount

On Order



[Ask a question about this product](#)

Manufacturer [Corning](#)

Reviews

There are yet no reviews for this product.