

10w/v% Polyoxyethylene(20) Sorbitan Monolaurate Solution

The Wako logo consists of the word "Wako" in white, sans-serif font, centered within a solid red square.The Wako logo consists of the word "Wako" in white, sans-serif font, centered within a solid red square.

100ml

Rating: Not Rated Yet

Price

Sales price 1080 ?

Discount

On Order

★★★★★

[Ask a question about this product](#)

Manufacturer [FUJIFILM Wako](#)

Reviews

There are yet no reviews for this product.