

## Wakopak Wakosil-II5C18RS Prep 4.6\*250mm

The Wako logo consists of the word "Wako" in white, sans-serif font, centered within a solid red square.The Wako logo consists of the word "Wako" in white, sans-serif font, centered within a solid red square.

1 column(D)

Rating: Not Rated Yet

**Price**

Sales price 4820 ?

Discount

*On Order*

★★★★★

[Ask a question about this product](#)

Manufacturer [FUJIFILM Wako](#)

### Reviews

There are yet no reviews for this product.