

## MagCaptur™ HP Anti PA tag Antibody Magnetic Beads

The Wako logo consists of the word "Wako" in white, sans-serif font, centered within a solid red square.The Wako logo consists of the word "Wako" in white, sans-serif font, centered within a solid red square.

2 ml

Rating: Not Rated Yet

**Price**

Sales price 7290 ?

Discount

[Ask a question about this product](#)

Manufacturer [FUJIFILM Wako](#)

Description

[MagCapture™ HP?PA tag????](#)

**Reviews**

There are yet no reviews for this product.